



External Carbon Dioxide Probe

Specification

1. Probe range: 0 to 5000 ppm
2. Probe accuracy: ± 50 ppm

What's in the Pack


CO₂ gas probe

Preparing the Probe

1. Connect the CO₂ probe adapter to the Micro-USB external probe socket on the Ward's DataHub. Make sure that the USB





sign on the adapter cable connector is facing up.

2. Turn on the Ward's DataHub.
3. Click on the probe key  located above the external probe socket where you've plugged in the CO₂ probe. If this key has two probe icons, you might need to press it twice to show the magnetic field probe.
4. The probe needs 30 seconds to "warm-up". During this time, the Ward's DataHub will show 5000 ppm. After the "warm-up" time the Ward's DataHub will display 350 to 600 ppm in a ventilated room.
5. Using the CO₂ probe, the Ward's DataHub will not enter "sleep-mode" to avoid warming up the probe for 30 seconds before each measurement.
6. The Ward's DataHub shows the CO₂ probe readings in "ppm" units.

Using the CO₂ Probe with the Ward's DataHub Analysis Software/App

1. Open the Ward's DataHub Analysis Software/App.
2. Make sure it is connected to the Ward's DataHub via USB or Bluetooth communication.



3. Click the SETUP icon  to open the *Logger setup* dialog box:
4. The CO₂ probe is automatically identified by the software.
5. To select this probe, simply click on selection button to the left of the CO₂ probe icon.