

# **Specification**

- 1. Probe range: -10 V to +10 V
- 2. Probe accuracy: ±2 %

### What's in the Pack

- ① Voltage probe adapter
- (2) A pair of red and black banana cables

#### Preparing the Probe

 Connect the Voltage probe adapter to the Micro-USB external probe socket on the Ward's DataHub. Make sure that the USB

sign V on the adapter cable connector is facing up.

Connect the red and black banana cables to the other end of the Voltage probe.

### <u>Using the Voltage Probe with the Ward's</u> <u>DataHub</u>

- 1. Turn on the Ward's DataHub.
- The Ward's DataHub reads the Voltage probe level. If it is within ±10 % around zero, the Ward's DataHub will perform an automatic probe offset calibration, showing a zero reading.
- 3. The other end of the banana cables can be connected to the experiment setup.

- 4. Click on the probe key located above the external probe socket where the Voltage probe is plugged in. If this key has two probe icons, it may need to be pressed twice to show the Voltage probe.
- The Ward's DataHub shows the voltage reading in voltage units.

## <u>Using the Voltage Probe with the Ward's</u> <u>DataHub Analysis Software</u>

- Open the Ward's DataHub Analysis Software.
- Make sure it is connected to the Ward's DataHub via USB or Bluetooth communication.



### 3. Click the SETUP icon

to open

Select Sensors	
Ar pressure	
🖻 🍯 Voltage	
🗉 🚎 Low Voltage	
E 🔆 Late	Rate
😢 🇭 Voltage	1/Sec 🔄
🗐 🔓 Amb. Temperature	Samples
🗉 🔐 Ext. Temperature	100 -
🗈 🍏 Current	
Distance	
🔲 뤷 W Microphone	Exit
🗈 💐 Acceleration	
🛅 🤯 Heart rate	
🗉 🌡 Barometer	

the Logger setup dialog box:

- The Voltage probe is automatically identified by the software.
- To select this probe, simply click on selection button to the left of the Voltage probe icon.