



EPISIOTOMY AND SUTURE TRAINER

For realistic training of professional care and suturing of birth injuries:

- Make individual episiotomy incisions and sutures
- Palpate muscle structures in the perineal area
- Treat vaginal and labial tears with usual suture techniques
- It is possible to suture an incision or tear multiple times



TRAINING ON INDIVIDUAL INCISIONS, SUCH AS:

- Median episiotomy
- Mediolateral episiotomy
- Lateral episiotomy

Thanks to the episiotomy and suture trainer, it is possible to demonstrate and practice professionally caring for and suturing birth injuries. The easy-to-clean material was specially developed for the training requirements of wound management. Incisions can be made and then sutured up again multiple times.



PALPATE THE FOLLOWING MUSCLE STRUCTURES IN THE PERINEAL AREA:

- M. sphincter ani externus
- M. transversus perinei superficialis
- M. bulbospongiosus



Muscle and skin layers in the perineal area that can be sutured independently



Muscle and skin layers can be sutured separately. The suture procedure for each incision can be carried out multiple times.

› Delivery contents:

Episiotomy and Suture Trainer

- Silicone insert for episiotomy training
- Base and stand for secure and angled positioning on a table

P95 / 1019639

› Order list:

Replacement set for episiotomy trainer

XP95-001 / 1019640

Replacement sets for episiotomy trainer (set of 5)

XP95-002 / 1019641

› Optional:

Suture set episiotomy and suture trainer

(sutures, needle holder, scissors, tweezers, scalpel holder, blades)

XP95-003 / 1020767

› Recommended:

Birth stages trainer with 6 stages

P94 / 1020628

Birth simulator PRO

P90 / 1020333

9000772 10/2016

All training products for midwifery and obstetrics can be found online at **3bscientific.com**



3B Scientific GmbH • Rudorffweg 8 • 21031 Hamburg • Germany

Tel.: +49 40 739 66-0 • **Fax:** +49 40 739 66-100 • **E-Mail:** info@3bscientific.com